VI Olimpíada de Biologia de Catalunya

Febrer 2016

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DNA Multi-Sample Kit

REAGENTS

- DNA Binding Beads, 10 mg/mL
- DNA Elution Buffer 1
- DNA Elution Buffer 2
- Multi-Sample DNA Lysis Buffer
- Wash Solution 1
- Wash Solution 2
- Water nuclease-free (DNA resuspension)
- Water sterile (mouth wash)

ADAPTED PROTOCOL FOR THE "OLIMPIADES DE BIOLOGIA 2016
(CATALUNYA)" BY M^a JOSÉ LÓPEZ-INIESTA AND GEMMA MARFANY
(UNIVERSITAT DE BARCELONA)

GENERAL CONSIDERATIONS

Guidelines

- Perform the protocol at room temperature (20 to 25 °C) except where noted.
- Store and use the DNA Elution Buffers at room temperature.
- Avoid creating bubbles when mixing samples by pipetting up/down.
- When aspirating, be careful not to dislodge the DNA binding beads from the magnet.
- When capturing beads on the magnetic stand, you can remove the supernatant once the solution is clear. The bead collection times can vary depending on sample type and the quantity of nucleic acid used.
- Close the tubes during the binding, washing, and elution steps to prevent spill-over and cross-contamination.
- The vortex shaking speeds in this protocol are recommendations, as the speeds can vary depending on the model of vortexer used. Ideal speeds mix the samples thoroughly without splashing.
- Use 2-mL reaction tubes for all samples.
- During the washing and elution steps, vortex the samples in short pulses to prevent the beads from sticking to the sides of the tubes.
- If the DNA Binding Beads are being removed with the supernatant, use smaller pipet tips to remove the supernatant

Isolation of genomic DNA from buccal swabs

Before you begin

- Review the guidelines described in the first page.
- Preheat a heated block (or alternate heat source) to 70 °C.
- Be sure you have all the reagents you need before starting the protocol.

Collect and disrupt the sample

1. Use 10 ml sterile water to wash thoroughly your mouth. Keep moving the water inside your buccal cavity back and forth for several minutes.
2. Collect at least 5 ml in a 15 ml plastic sterile tube.
3. Centrifuge at 2000 rpm for 5 min. The buccal mucose cells will pellet at the bottom.
4. Discard the supernatant and resuspend with For each buccal swab, add 200 μ L of Multi-Sample DNA Lysis Buffer (MS tube)
5. Collect the cells and pass them into a fresh 2 mL processing tube.
6. Shake the tube for 3 minutes at speed 1 or 2 on the vortex adaptor, until the pellet is resuspended.

Perform the DNA extraction

1. Add 160 μ L of 100% isopropanol to each lysate (IP tube)
2. Shake end-over-end the processing tube vigorously for 3 minutes.
3. Add 40 μ L of prepared DNA-binding bead mix to the sample, then shake the tube end-over-end with vigour (Mag tube)
4. Place the tube on the magnetic stand for 2 minutes or until the solution clears and the beads are pelleted against the magnet.
5. Discard the supernatant carefully, avoiding the beads.
6. Wash the sample using 500 μ L of Wash Solution 1 (WS1 tube). Take the tube out of the magnet and shake vigorously end-over-end.
7. Place the tube on the magnetic stand for 1-2 minutes until the beads are pelleted against the magnet.
8. Discard the supernatant being careful not to disturb the beads. Repeat step 5 using 500 μ L of Wash Solution 2 (WS2 tube)
9. Keeping the tube on the magnetic stand, discard the supernatant, being careful not to disturb the beads. Air-dry the uncapped tubes on the magnetic stand for 3 minutes.

10. OPTIONAL!!!! OTHERWISE SKIP AND GO TO POINT 11

- While the tubes are drying, prepare sufficient RNase A mix for your samples.
- Remove the tube from the magnetic stand, add 100 μ L of RNase mixture to the sample, then shake the tube for 2 minutes at speed 1 or 2 on the vortex adaptor.
- Remove the tube from the vortex adaptor, add 100 μ L of Multi-Sample DNA Lysis Buffer and 120 μ L of 100% isopropanol to the sample, then shake the tube for 3 minutes at the vortex.
- Place the tube on the magnetic stand for 1 minute until the beads are pelleted against the magnet.
- Repeat step 5 twice more using Wash Solution 2. While the tube is on the magnet, discard the supernatant, then air-dry the uncapped sample on the magnetic stand for 3 minutes

Perform the elution

11. Add 25 μ L of DNA Elution Buffer 1 to the sample (EB1 tube)
12. Briefly vortex the sample to break apart the beads.
13. Incubate the sample at 70 °C for 5 min on a heated block.
14. Remove the tube and add 25 μ L of DNA Elution Buffer 2 to the sample (EB2 tube)
15. Gently vortex the sample for 20 to 30 seconds until resuspended, then place the tube on the magnetic stand for 5 minutes or until the solution clears and the beads are pelleted against the magnet.
16. While the tube is on the magnet, transfer the eluate (which contains the purified DNA) to a new tube, then close the tube immediately.

DNA Electrophoresis and Visualization

17. Take 16 μ L of your DNA and add 4 μ L of Loading Buffer x 5 (TC tube).
18. Load it into the well of an horizontal agarose gel.
19. Connect the electrophoresis source to 100 V. Let it run for 30 min
20. Visualize it under the UV light.

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